



pth Antibody

Product Code	CSB-PA358958HA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7D1
Immunogen	Recombinant Escherichia coli Peptidyl-tRNA hydrolase protein (1-194AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-tRNA hydrolase (PTH) (EC 3.1.1.29), pth
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	pth
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.