



recF Antibody, HRP conjugated

Product Code	CSB-PA358973LB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7H0
Immunogen	Recombinant Escherichia coli DNA replication and repair protein RecF protein (2-357AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA replication and repair protein RecF, recF, uvrF
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	recF