



recF Antibody, Biotin conjugated

Product Code	CSB-PA358973LD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7H0
Immunogen	Recombinant Escherichia coli DNA replication and repair protein RecF protein (2-357AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA replication and repair protein RecF, recF, uvrF
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	recF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.