



ruvA Antibody, FITC conjugated

Product Code	CSB-PA359062HC01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A809
Immunogen	Recombinant Escherichia coli Holliday junction ATP-dependent DNA helicase RuvA protein (1-203AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Holliday junction ATP-dependent DNA helicase RuvA (EC 3.6.4.12), ruvA
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	ruvA