



fabF Antibody, FITC conjugated

Product Code	CSB-PA359466LC01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AAI5
Immunogen	Recombinant Escherichia coli 3-oxoacyl-[acyl-carrier-protein] synthase 2 protein (2-413AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-oxoacyl-[acyl-carrier-protein] synthase 2 (EC 2.3.1.179) (3-oxoacyl-[acyl-carrier-protein] synthase II) (Beta-ketoacyl-ACP synthase II) (KAS II), fabF, fabJ
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	fabF