



fabF Antibody, Biotin conjugated

Product Code	CSB-PA359466LD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AAI5
Immunogen	Recombinant Escherichia coli 3-oxoacyl-[acyl-carrier-protein] synthase 2 protein (2-413AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	3-oxoacyl-[acyl-carrier-protein] synthase 2 (EC 2.3.1.179) (3-oxoacyl-[acyl-carrier-protein] synthase II) (Beta-ketoacyl-ACP synthase II) (KAS II), fabF, fabJ
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	fabF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.