



repD Antibody, FITC conjugated

Product Code	CSB-PA360946HC01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03065
Immunogen	Recombinant Staphylococcus aureus Replication initiation protein (1-311AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Replication initiation protein, repD
Immunogen Species	Staphylococcus aureus
Research Area	Others
Target Names	repD