



repD Antibody, Biotin conjugated

Product Code	CSB-PA360946HD01FKZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03065
Immunogen	Recombinant Staphylococcus aureus Replication initiation protein (1-311AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Replication initiation protein, repD
Immunogen Species	Staphylococcus aureus
Research Area	Others
Target Names	repD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.