



BLLF1 Antibody, Biotin conjugated

Product Code	CSB-PA360986LD01EFA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03200
Immunogen	Recombinant Epstein-Barr virus Envelope glycoprotein GP350 protein (1-493AA)
Raised In	Rabbit
Species Reactivity	Epstein-Barr virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein GP350 (Membrane antigen) (MA)
Immunogen Species	Epstein-Barr virus (strain B95-8) (HHV-4) (Human herpesvirus 4)
Research Area	Others
Target Names	BLLF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.