



Trypsin inhibitor DE-3 Antibody, Biotin conjugated

Product Code	CSB-PA362840YD01EGO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09943
Immunogen	Recombinant Erythrina caffra Trypsin inhibitor DE-3 protein (1-172AA)
Raised In	Rabbit
Species Reactivity	Erythrina caffra
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trypsin inhibitor DE-3
Immunogen Species	Erythrina caffra (Kaffir tree) (Coastal coral tree)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.