



gshA Antibody, Biotin conjugated

Product Code	CSB-PA363915HD01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A6W9
Immunogen	Recombinant Escherichia coli Glutamate--cysteine ligase protein (1-518AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamate--cysteine ligase (EC 6.3.2.2) (Gamma-ECS) (GCS) (Gamma-glutamylcysteine synthetase), gshA, gsh-I
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	gshA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.