



# recO Antibody, HRP conjugated

Product Code	CSB-PA364006HB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A7H3
Immunogen	Recombinant Escherichia coli DNA repair protein RecO protein (1-242AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA repair protein RecO (Recombination protein O), recO
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	recO