



recO Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA364006HC01ENV |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P0A7H3 |
| Immunogen | Recombinant Escherichia coli DNA repair protein RecO protein (1-242AA) |
| Raised In | Rabbit |
| Species Reactivity | Escherichia coli |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | DNA repair protein RecO (Recombination protein O), recO |
| Immunogen Species | Escherichia coli (strain K12) |
| Research Area | Others |
| Target Names | recO |