



ruvC Antibody

Product Code	CSB-PA364095HA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A814
Immunogen	Recombinant Escherichia coli Crossover junction endodeoxyribonuclease RuvC protein (2-173AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Crossover junction endodeoxyribonuclease RuvC (EC 3.1.22.4) (Holliday junction nuclease RuvC) (Holliday junction resolvase RuvC), ruvC
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	ruvC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.