



# sulA Antibody, HRP conjugated

Product Code	CSB-PA365383LB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0AFZ5
Immunogen	Recombinant Escherichia coli Cell division inhibitor SulA protein (1-169AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division inhibitor SulA, sulA, sfiA
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	sulA