



uvsY Antibody, HRP conjugated

Product Code	CSB-PA366179LB01EDZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04537
Immunogen	Recombinant Enterobacteria phage T4 Recombination protein uvsY protein (1-137AA)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage T4
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Recombination protein uvsY, uvsY
Immunogen Species	Enterobacteria phage T4 (Bacteriophage T4)
Research Area	Others
Target Names	uvsY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.