



PR10A Antibody, FITC conjugated

Product Code	CSB-PA381182LC01CXF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2A1A1
Immunogen	Recombinant Coptis japonica S-norcochlorine synthase 2 protein (20-196AA)
Raised In	Rabbit
Species Reactivity	Coptis japonica
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	S-norcochlorine synthase 2 (EC 4.2.1.78) (Pathogenesis related protein 10A) (CjPR10A), PR10A
Immunogen Species	Coptis japonica (Japanese goldthread)
Research Area	Others
Target Names	PR10A