



# amiF Antibody

|                     |  |
|---------------------|--|
| Product Code        | CSB-PA397462LA01BPF  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | A4Z3G9   |
| Immunogen           | Recombinant Bradyrhizobium sp. Formamidase protein (1-337AA)                     |
| Raised In           | Rabbit   |
| Species Reactivity  | Bradyrhizobium sp.   |
| Tested Applications | ELISA  |
| Form                | Liquid   |
| Conjugate           | Non-conjugated   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified   |
| Isotype             | IgG  |
| Clonality           | Polyclonal   |
| Alias               | Formamidase (EC 3.5.1.49) (Formamide amidohydrolase), amiF                       |
| Immunogen Species   | Bradyrhizobium sp. (strain ORS 278)  |
| Research Area       | Others   |
| Target Names        | amiF   |