



amiF Antibody, HRP conjugated

Product Code	CSB-PA397462LB01BPF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4Z3G9
Immunogen	Recombinant Bradyrhizobium sp. Formamidase protein (1-337AA)
Raised In	Rabbit
Species Reactivity	Bradyrhizobium sp.
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Formamidase (EC 3.5.1.49) (Formamide amidohydrolase), amiF
Immunogen Species	Bradyrhizobium sp. (strain ORS 278)
Research Area	Others
Target Names	amiF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.