



amiF Antibody, Biotin conjugated

Product Code	CSB-PA397462LD01BPF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4Z3G9
Immunogen	Recombinant Bradyrhizobium sp. Formamidase protein (1-337AA)
Raised In	Rabbit
Species Reactivity	Bradyrhizobium sp.
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Formamidase (EC 3.5.1.49) (Formamide amidohydrolase), amiF
Immunogen Species	Bradyrhizobium sp. (strain ORS 278)
Research Area	Others
Target Names	amiF