



mip Antibody, Biotin conjugated

Product Code	CSB-PA401660HD01LLJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A5IGB8
Immunogen	Recombinant Legionella pneumophila Outer membrane protein MIP protein (21-233AA)
Raised In	Rabbit
Species Reactivity	Legionella pneumophila
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein MIP (EC 5.2.1.8) (Macrophage infectivity potentiator) (Peptidyl-prolyl cis-trans isomerase) (PPIase) (Rotamase), mip
Immunogen Species	Legionella pneumophila (strain Corby)
Research Area	Others
Target Names	mip
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.