



# egsA Antibody, FITC conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA455307HC01CFM  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | B3EYN8   |
| <b>Immunogen</b>           | Recombinant Geobacillus stearothermophilus Glycerol-1-phosphate dehydrogenase [NAD(P)+] protein (1-404AA)  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Geobacillus stearothermophilus   |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | FITC   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4   |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Glycerol-1-phosphate dehydrogenase [NAD(P)+] (G1P dehydrogenase) (G1PDH) (EC 1.1.1.261) (Enantiomeric glycerophosphate synthase) (sn-glycerol-1-phosphate dehydrogenase), egsA, araM |
| <b>Immunogen Species</b>   | Geobacillus stearothermophilus (Bacillus stearothermophilus)   |
| <b>Research Area</b>       | Others   |
| <b>Target Names</b>        | egsA   |