



# egsA Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA455307HD01CFM
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	B3EYN8
<b>Immunogen</b>	Recombinant Geobacillus stearothermophilus Glycerol-1-phosphate dehydrogenase [NAD(P)+] protein (1-404AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Geobacillus stearothermophilus
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Glycerol-1-phosphate dehydrogenase [NAD(P)+] (G1P dehydrogenase) (G1PDH) (EC 1.1.1.261) (Enantiomeric glycerophosphate synthase) (sn-glycerol-1-phosphate dehydrogenase), egsA, araM
<b>Immunogen Species</b>	Geobacillus stearothermophilus (Bacillus stearothermophilus)
<b>Research Area</b>	Others
<b>Target Names</b>	egsA