



tig Antibody, Biotin conjugated

Product Code	CSB-PA484773LD01EOB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7UJQ9
Immunogen	Recombinant Escherichia coli O127:H6 Trigger factor protein (1-432AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O127:H6
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Trigger factor (TF) (EC 5.2.1.8) (PPlase), tig
Immunogen Species	Escherichia coli O127:H6 (strain E2348/69 / EPEC)
Research Area	Others
Target Names	tig
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.