



blaNDM-1 Antibody, Biotin conjugated

Product Code	CSB-PA510513HD01KBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C7C422
Immunogen	Recombinant Klebsiella pneumoniae Metallo-beta-lactamase type 2 protein (29-270AA)
Raised In	Rabbit
Species Reactivity	Klebsiella pneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metallo-beta-lactamase type 2 (EC 3.5.2.6) (B2 metallo-beta-lactamase) (Beta-lactamase type II) (Metallo-beta-lactamase NDM-1) (Metallo-beta-lactamase type II) (New Delhi metallo-beta-lactamase-1) (NDM-1), blaNDM-1
Immunogen Species	Klebsiella pneumoniae
Research Area	Others
Target Names	blaNDM-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.