



aiiM Antibody, Biotin conjugated

Product Code	CSB-PA512020HD01MTN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C6L862
Immunogen	Recombinant Microbacterium testaceum N-acyl homoserine lactonase protein (1-251AA)
Raised In	Rabbit
Species Reactivity	Microbacterium testaceum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acyl homoserine lactonase (AHL-lactonase) (EC 3.1.1.81), aiiM
Immunogen Species	Microbacterium testaceum (strain StLB037)
Research Area	Others
Target Names	aiiM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.