



RPN13 Antibody, Biotin conjugated

Product Code	CSB-PA517092HD01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O13563
Immunogen	Recombinant Saccharomyces cerevisiae 26S proteasome regulatory subunit RPN13 protein (2-156AA)
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	26S proteasome regulatory subunit RPN13 (Proteasome non-ATPase subunit 13), RPN13, DAQ1
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Research Area	Others
Target Names	RPN13
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.