



IspL Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA526307HA01RKU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | O54067 |
| Immunogen | Recombinant Rhizobium meliloti UDP-glucuronate 5'-epimerase protein (1-341AA) |
| Raised In | Rabbit |
| Species Reactivity | Rhizobium meliloti |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | UDP-glucuronate 5'-epimerase (EC 5.1.3.12) (UDP-glucuronic acid epimerase), IspL |
| Immunogen Species | Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti) |
| Research Area | Others |
| Target Names | IspL |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |