



IspL Antibody, Biotin conjugated

Product Code	CSB-PA526307HD01RKU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O54067
Immunogen	Recombinant Rhizobium meliloti UDP-glucuronate 5'-epimerase protein (1-341AA)
Raised In	Rabbit
Species Reactivity	Rhizobium meliloti
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-glucuronate 5'-epimerase (EC 5.1.3.12) (UDP-glucuronic acid epimerase), IspL
Immunogen Species	Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti)
Research Area	Others
Target Names	IspL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.