



lepB Antibody, HRP conjugated

Product Code	CSB-PA604637LB01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WKA1
Immunogen	Recombinant Mycobacterium tuberculosis Signal peptidase I protein (88-294AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptidase I (SPase I) (EC 3.4.21.89) (Leader peptidase I), lepB
Immunogen Species	Mycobacterium tuberculosis
Research Area	Others
Target Names	lepB