



lepB Antibody, HRP conjugated

Product Code	CSB-PA604637LB01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WKA1
Immunogen	Recombinant Mycobacterium tuberculosis Signal peptidase I protein (88-294AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptidase I (SPase I) (EC 3.4.21.89) (Leader peptidase I), lepB
Immunogen Species	Mycobacterium tuberculosis
Research Area	Others
Target Names	lepB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.