




INCA1 Antibody

Product Code	CSB-PA605764LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0VD86
Immunogen	Recombinant Human Protein INCA1 protein (1-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein INCA1 (Inhibitor of CDK interacting with cyclin A1), INCA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	INCA1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA605764LA01HU at dilution of 1:100