



INCA1 Antibody, Biotin conjugated

Product Code	CSB-PA605764LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0VD86
Immunogen	Recombinant Human Protein INCA1 protein (1-236AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein INCA1 (Inhibitor of CDK interacting with cyclin A1), INCA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	INCA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.