



RASGEF1B Antibody, Biotin conjugated

Product Code	CSB-PA610807LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0VAM2
Immunogen	Recombinant Human Ras-GEF domain-containing family member 1B protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-GEF domain-containing family member 1B (GPI gamma-4), RASGEF1B, GPIG4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASGEF1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.