



TMEM267 Antibody, HRP conjugated

Product Code	CSB-PA610832LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0VDI3
Immunogen	Recombinant Human Transmembrane protein 267 protein (1-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 267, TMEM267, C5orf28
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM267
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.