



IL24 Antibody, Biotin conjugated

Product Code	CSB-PA613378LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13007
Immunogen	Recombinant Human Interleukin-24 protein (52-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-24 (IL-24) (Melanoma differentiation-associated gene 7 protein) (MDA-7) (Suppression of tumorigenicity 16 protein), IL24, MDA7 ST16
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.