



SEMA5A Antibody, HRP conjugated

Product Code	CSB-PA613413LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13591
Immunogen	Recombinant Human Semaphorin-5A protein (201-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Semaphorin-5A (Semaphorin-F) (Sema F), SEMA5A, SEMAF
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SEMA5A