



LY6D Antibody

Product Code	CSB-PA613492LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14210
Immunogen	Recombinant Human Lymphocyte antigen 6D protein (21-98AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte antigen 6D (Ly-6D) (E48 antigen), LY6D, E48
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LY6D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.