



KCNMA1 Antibody, FITC conjugated

Product Code	CSB-PA614255LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12791
Immunogen	Recombinant Human Calcium-activated potassium channel subunit alpha-1 protein (1-86AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-activated potassium channel subunit alpha-1 (BK channel) (BKCA alpha) (Calcium-activated potassium channel, subfamily M subunit alpha-1) (K(VCA)alpha) (KCa1.1) (Maxi K channel) (MaxiK) (Slo-alpha) (Slo1) (Slowpoke homolog) (Slo homolog) (hSlo), KCNMA1, KCNMA SLO
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	KCNMA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.