



GPM6B Antibody, Biotin conjugated

Product Code	CSB-PA614417LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13491
Immunogen	Recombinant Human Neuronal membrane glycoprotein M6-b protein (1-64AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neuronal membrane glycoprotein M6-b (M6b), GPM6B, M6B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPM6B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.