



POLA2 Antibody, HRP conjugated

Product Code	CSB-PA614517LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14181
Immunogen	Recombinant Human DNA polymerase alpha subunit B protein (101-400AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase alpha subunit B (DNA polymerase alpha 70 kDa subunit), POLA2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.