



ZMYM3 Antibody, HRP conjugated

Product Code	CSB-PA614520LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14202
Immunogen	Recombinant Human Zinc finger MYM-type protein 3 protein (1-495AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger MYM-type protein 3 (Zinc finger protein 261), ZMYM3, DXS6673E KIAA0385 ZNF261
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZMYM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.