



UBXN2B Antibody, Biotin conjugated

Product Code	CSB-PA614827LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14CS0
Immunogen	Recombinant Human UBX domain-containing protein 2B protein (2-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UBX domain-containing protein 2B (NSFL1 cofactor p37) (p97 cofactor p37), UBXN2B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UBXN2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.