



LONRF1 Antibody, Biotin conjugated

Product Code	CSB-PA615067LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q17RB8
Immunogen	Recombinant Human LON peptidase N-terminal domain and RING finger protein 1 protein (538-773AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LON peptidase N-terminal domain and RING finger protein 1 (RING finger protein 191), LONRF1, RNF191
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LONRF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.