



SMT3 Antibody, HRP conjugated

Product Code	CSB-PA615489LB01SVG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q12306
Immunogen	Recombinant Saccharomyces cerevisiae Ubiquitin-like protein SMT3 protein (2-98AA)
Raised In	Rabbit
Species Reactivity	Saccharomyces cerevisiae
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-like protein SMT3, SMT3
Immunogen Species	Saccharomyces cerevisiae (strain ATCC 204508 / S288c) (Baker's yeast)
Research Area	Cell Biology
Target Names	SMT3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.