



PIGH Antibody, HRP conjugated

Product Code	CSB-PA617916LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14442
Immunogen	Recombinant Human Phosphatidylinositol N-acetylglucosaminyltransferase subunit H protein (1-188AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol N-acetylglucosaminyltransferase subunit H (EC 2.4.1.198) (Phosphatidylinositol-glycan biosynthesis class H protein) (PIG-H), PIGH
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PIGH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.