



NFIX Antibody, Biotin conjugated

Product Code	CSB-PA617941LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14938
Immunogen	Recombinant Human Nuclear factor 1 X-type protein (201-440AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear factor 1 X-type (NF1-X) (Nuclear factor 1/X) (CCAAT-box-binding transcription factor) (CTF) (Nuclear factor I/X) (NF-I/X) (NFI-X) (TGGCA-binding protein), NFIX
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NFIX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.