



SPCS2 Antibody, Biotin conjugated

Product Code	CSB-PA617993LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15005
Immunogen	Recombinant Human Signal peptidase complex subunit 2 protein (7-86AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptidase complex subunit 2 (EC 3.4.-.-) (Microsomal signal peptidase 25 kDa subunit) (SPase 25 kDa subunit), SPCS2, KIAA0102 SPC25
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SPCS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.