



SEPT2 Antibody, Biotin conjugated

Product Code	CSB-PA617994LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15019
Immunogen	Recombinant Human Septin-2 protein (1-361AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Septin-2 (Neural precursor cell expressed developmentally down-regulated protein 5) (NEDD-5), SEPT2, DIFF6 KIAA0158 NEDD5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SEPT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.