



EPHA7 Antibody, FITC conjugated

Product Code	CSB-PA618004LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15375
Immunogen	Recombinant Human Ephrin type-A receptor 7 protein (28-182AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ephrin type-A receptor 7 (EC 2.7.10.1) (EPH homology kinase 3) (EHK-3) (EPH-like kinase 11) (EK11) (hEK11), EPHA7, EHK3 HEK11
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EPHA7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.